


# Secondary Ion Mass Spectrometry

Ion TOF.SIMS<sup>5</sup>

## PRINCIPLE

Time-of-Flight Secondary Ion Mass Spectrometry (TOF-SIMS) is a very sensitive surface analytical technique which provides detailed elemental information about the surface, thin layers, interfaces of the sample, and gives a full three-dimensional analysis. Finely focused ion beam sputters sample surface and the exact mass of emitted ions and ion clusters is measured by time-of-flight analyzer. From the exact mass and intensity of the SIMS peak, the identity of an element or molecular fragments can be determined.


## SCHEMA



## SPECIFICATION

<b>Primary ions</b>	Bi <sup>1+</sup> , Bi <sup>1++</sup> , Bi <sup>2+</sup> .. Bi <sub>z</sub> <sup>+</sup> , Mn <sup>+</sup> Energy: 30 keV
<b>Sputtering ions</b>	Cs <sup>+</sup> , O <sub>2</sub> <sup>+</sup> Energy: 0.5 – 2.0 keV
<b>Electron Flood gun</b>	Energy: <20 eV
<b>Sample holders</b>	Back side mounting stage: Sample size: 15 mm × 10 mm Top side mounting stage: Sample size up to 100 mm × 50 mm Heating/Cooling stage: Sample size: 10 mm × 10 mm Temp. range -130°C – 600°C Rotating stage: Sample size: Ø10 mm
<b>Surface spectrometer</b>	High sensitivity (1–2 ML)
<b>Surface imaging</b>	High lateral resolution (<60 nm)
<b>Depth profiling</b>	Depth resolution better than 1 nm High mass resolution > 11 000 @ 29 u Sputter speed of up to 10 µm/h
<b>3D Analysis</b>	Parallel mass detection

## RESULTS



left – High resolution mass spectrum of a PET sample.  
right – Mass resolved 3D analysis of TFT display pixel.

## MORE INFO

**Guarantor:** Marek Otevřel (marek.otevrel@ceitec.vutbr.cz)

**Web:** <http://nano.ceitec.cz/secondary-ion-mass-spectroscopy-ion-tof-tof-sims5-sims/>

